



LCO₂ Filter Housing (# 8784) and High Pressure Hose (# 8770-32, # 8770-33, # 8770-HPS)



8784 LCO₂ Filter Housing with # 8760-34 CGA - 320 S.S. Coupler & # 8770-33 High Pressure Hose (10ft)

- ① Gasket (Teflon or Nylon) (# 8784-05)
- ② CGA - 320 S.S. Coupler Connection Nut to LCO₂ Tank (# 8760-34: ② + ③)
- ③ CGA - 320 S.S. Nipple (# 8760-34: ② + ③)
- ④ Filter Housing w/ Oil/Water Filter Element (# 8784A)
- ⑤ T-Filter w/ 0.5 µm Particulate Element (# 8770-83B)
- ⑥ High Pressure Hose to Samdri® (Sold separately: 5ft - # 8770-32, 10ft - # 8770-33, Custom - # 8770-HPS)



8784 LCO₂ Filter Housing (④ + ⑤)

High Pressure Hose to Samdri®



Connector to SAMDRI®