


8785 Double T-Filter Assembly



1 Main LCO₂ Tank Valve

2 Teflon® Gasket (# 8784-05)

3 CGA-320 S.S. Nut (# 8760-34: ③ + ⑥)

4 T-Filter Assembly #2 with T-Filter Element (0.5 µm, # 8770-83B)

5 T-Filter Assembly #1 with T-Filter Element (0.5 µm, #8770-83B)

6 CGA-320 S.S. Coupler Nipple (# 8760-34: ③ + ⑥)

7 Teflon® tape

8 High Pressure Hose to Samdri® (Sold separately: 5ft - # 8770-32, 10ft - # 8770-33, Custom - # 8770-HPS)